



tebubio

Facilitators of Life Sciences Research

# How to Identify Efficient LNP Formulations for Your Model ?

## Insights from Tebubio's Delivery Platform

Xavier Warnet, PhD – Project Manager

February 6<sup>th</sup> 2025



# Our Team Today



## Speaker

Xavier Warnet, PhD

Project Manager, Formulation Specialist



## Q&A Moderator

Erica Cirri, PhD

Project Manager, RNA Specialist



## Session Moderator

Frédéric Samazan

Event Manager

# Webinar Agenda

1. | Tebubio at a glance
2. | Lipid NanoParticles basics
3. | Delivery & Formulation Challenges
4. | Overview of Tebubio Delivery Platform
5. | Case Study: LNP formulation screening for Cancer Cells Delivery
6. | Innovations through our parnters
7. | Live Q&A Session

# Webinar Agenda

1. Tebubio at a glance
2. Lipid NanoParticles basics
3. Delivery & Formulation Challenges
4. Overview of Tebubio Delivery Platform
5. Case Study: LNP formulation screening for Cancer Cells Delivery
6. Innovations through our parnters
7. Live Q&A Session

# We facilitate Life Sciences Research everyday

and contribute to a brighter future



## We are Pan-European

- Founded in 1953
- Family-owned
- 100+ Employees
- Local offices across Europe



## We Act for Life Sciences

- Innovation is in our DNA
- Contract Research Services Lab
- Part of EU Life Sciences ecosystems



## We Care

- Ethical, compliant and transparent sourcing (from OEMs only)
- Animal welfare policy
- Corporate Social Responsibility (ISO 14001 / ISO 9001, Decarbonation Program)

# With Tebubio, advance your Research Projects faster

Thanks to our Holistic Range of Solutions



Order advanced biological solutions



Outsource and Accelerate research



Streamline your ordering process

## A large portfolio

Access to **over 1,300,000 standard** references and **non-catalogue** ones.

## From trusted & ethical suppliers

Get solutions from **reputable global suppliers**, carefully selected for their **relevance, ethical** and **legal** compliance (e.g. Animal Welfare).

## Dedicated scientific support

Our **Scientific Team** guides you to source, select and use solutions.

## Based in Europe

Our **Teams** and **Contract Research Services Lab** are in Europe.

## Team committed to success

A **PhD project manager** ensures the success of your project from A to Z.

## Strong expertise in Life sciences

- **mRNA** production & delivery
- **Cell line** engineering & **protein** production
- **Cellular studies**
- **Biomarkers** & Biostatistics analysis

## Reliable Supply Chain Management

- **IATA**
- **Human/Animal Biological Solutions**
- **Sourcing outside** existing suppliers
- **Warehousing** services

## Order from a single source

**Consolidate** your orders with us.

## Tailored agreements

From specific **one-off terms** to **procure-to-pay**, supported by **e-procurement** solutions.

# Tebubio CRS : Facilitators of Life Sciences Research

What do you want to talk about?



RNA based therapeutic  
discovery

Production | Delivery | Expression



*in vitro* modeling

Organoids | Microfluidics | 2&3D  
Models



Biomarker mapping

Target identification | Profiling |  
Quantification



Data analysis

Biostatistics | NGS | Custom

# Our Speaker



Xavier Warnet, PhD

CRS Laboratory Project Manager

- Xavier joined Teubio in February 2024 as a Project Manager on the RNA based Therapeutic Discovery platform . With a robust academic background, including a Ph.D. in Biochemistry and Biophysics from University Paris Cité, he brings deep expertise in biomolecular research, analytical techniques, and innovative drug delivery technologies.
- Prior to his role at Teubio, Xavier led the diagnostic team at LPS-BioSciences, focusing on the structural and functional characterization of Lipopolysaccharides (LPS) through advanced analytical techniques such as mass spectrometry (LC-MS-MS), LAL, and MAT assays.
- At Teubio, Xavier now leverages his expertise to guide clients in optimizing their delivery projects, and continues to develop and expand Teubio's LNP (Lipid Nanoparticle) delivery platform.



# Webinar Agenda





1. | Tebubio at a glance
2. | **Lipid NanoParticles basics**
3. | Delivery & Formulation Challenges
4. | Overview of Tebubio Delivery Platform
5. | Case Study: LNP formulation screening for Cancer Cells Delivery
6. | Innovations through our parnters
7. | Live Q&A Session


# Delivery or Formulation Challenges

## Some Basics on LNPs

- **Approved by the FDA for mRNA delivery.**
- **Biocompatible by nature.**
- **“Easy” to synthesize.**
- **Encapsulation of various cargos**  
(miRNA, siRNA, mRNA, ...)

LNPs are composed of (at least) 4 different "classes" of lipids:

- **Ionizable/cationic lipid (~50%):** 
- **Sterol (~38%):** 
- **Structural lipid (~10 %):** 
- **PEG-lipid (~1-2%):** 



Lipid nanoparticle (LNP)


*Hald Albertsen C. et al., Advanced Drug Delivery Reviews (2022)*

# Webinar Agenda

1. | Tebubio at a glance
2. | Lipid NanoParticles basics
3. | **Delivery & Formulation Challenges**
4. | Overview of Tebubio Delivery Platform
5. | Case Study: LNP formulation screening for Cancer Cells Delivery
6. | Innovations through our parnters
7. | Live Q&A Session

# Main obstacles for efficient LNP preparation and Delivery

## Delivery Services




**Formulation screening remains essential**

# Webinar Agenda

1. | Tebubio at a glance
2. | Lipid NanoParticles basics
3. | Delivery & Formulation Challenges
4. | **Overview of Tebubio Delivery Platform**
5. | Case Study: LNP formulation screening for Cancer Cells Delivery
6. | Innovations through our parnters
7. | Live Q&A Session

# We can Support You Throughout Your Entire mRNA Workflow

More insights on our Delivery Platform



# We can Support You Throughout Your Entire mRNA Workflow

What about custom mRNA production.



Want more information on our mRNA production platform and capabilities ?

Ask to receive the recording of our last webinar



- screening, to larger batch encapsulation for in vivo studies.
- **Accessibility** to a wide variety of off-the-shelf formulation and lipids.

# A complete solution for preclinical research

Our microfluidic devices

## For screening and initial optimization ...

Minimum amounts of reagents

Encapsulation from 10 $\mu$ g



*Spark NanoAssembler*  
Precision NanoSystems, Cytiva

## ... To advanced optimization and validation

Fine tuning of particles properties

Encapsulation to few mg RNA



*Tamara*  
Inside Therapeutics



# Webinar Agenda

1. | Tebubio at a glance
2. | Lipid NanoParticles basics
3. | Delivery & Formulation Challenges
4. | Overview of Tebubio Delivery Platform
5. | **Case Study: LNP formulation screening for Cancer Cells Delivery**
6. | Innovations through our parnters
7. | Live Q&A Session

# Case Study: LNP formulation screening for Cancer Cells Delivery

## Project


Identify the optimal LNP-based formulation to ensure efficient delivery of the candidate RNA to a spheroid cancer cell model (HCT116).

## Solution

- **15 different formulations** were settled, thanks to our large library of lipids, to identify the most efficient combination
- **Small scale screening** were performed, to identify a primary **selection of potent formulation specific** to the customer cellular models, including:
  - Quality control of LNPs physicochemical parameters : Encapsulation efficiency, average particle size and polydispersity index
  - Control of LNP Internalization and RNA translation efficiency : Fluorescent lipids and eGFP-mRNA.

## Customer pain point

- Finding the adapted lipids and formulation parameters.
- Characterize, select and optimize the best formulation.
- Assess transfection efficiency




# Case Study: LNP formulation screening for Cancer Cells Delivery

## LNP Formulation Design and Encapsulation



### 15 different formulations built

- Combining the most relevant lipids for the HCT 116 model
- Selection from our library and partners
- Formulations composed by:
  - **Mix of 5 different ionizable/cationic lipids**  
(ALC-0315, SM-102, LipidBrick IM21.7c, CP-LC-0741, DODMA)
  - **Mix of different structural lipids**  
(DSPC, DOPE)
  - **7 formulations with custom fluorescent lipids**  
(Cy3-1,2-Dipalmitoyl-sn-glycero-3-phosphoethanolamine)



### Encapsulation

- 5-15 µg of mRNA used for each formulation
- eGFP mRNA produced and optimised by Tebubio



# Case Study: LNP formulation screening for Cancer Cells Delivery

Quality control of LNPs physicochemical parameters




## Quality control

- **Encapsulation efficiency (EE)** analysis (%) by Ribogreen Assay
- **Average particle size (z-average)** by DLS
- **Polydispersity index (P)** by DLS

	1	2 *	3 *	4	5	6	7 *	8
EE (%)	98,2	62,5	41,3	63,3	99,6	100,0	43,4	87,7
z-average (nm)	84,4	75,0	72,1	74,8	120,1	109,8	77,4	65,6
PI	0,23	0,21	0,15	0,14	0,23	0,20	0,15	0,17
	9 *	10 *	11	12	13	14	15	
EE (%)	99,6	99,8	90,3	84,6	82,9	87,2	91,1	
z-average (nm)	95,4	100,4	90,7	104,3	129,9	100,1	95,7	
PI	0,29	0,24	0,20	0,07	0,06	0,07	0,09	

Formulations incorporating **fluorescent lipids**.



## Analysis of physicochemical parameters provide already some insights

- Poor encapsulation efficiency (\*): #3 #7 #2 & #4 → If those formulations proved to be efficient, more optimization would be required to achieve proper encapsulation.
- High heterogeneity (\*): #9 & #10 → which would be a problem for *in vivo* transfections but should have less of an impact for *in vitro* studies.


# Case Study: LNP formulation screening for Cancer Cells Delivery

Control of internalization and RNA expression



## Quality Control

- Transfection into HCT116 2D model
- 200ng of RNA used
- Analysis by microscopy : Brightfield (BF) (cell viability) and Fluorescence (eGFP) (translation level)



## Expression analysis allow us to perform a first selection of efficient formulations

- All the formulations tested led to protein expression.
- From this first screening, it was possible to select the most promising formulations : #4 #5 #7 & #8).
- It should be noted that poor encapsulation efficiencies is not necessarily connected to low protein expression (i.e #7, 43% encapsulation).


# Case Study: LNP formulation screening for Cancer Cells Delivery

Control of internalization and RNA expression



## Quality Control

- Transfection into HCT116 2D model
- 200ng of RNA used
- **Analysis by microscopy** :  
Brightfield (BF) (cell viability) and  
Fluorescence (translation level)



## Expression analysis allow us to perform a first selection of efficient formulations

- Fluorescent LNPs offer a complementary analysis.
- Membrane fusion **appears** to work for all of them (Cy3) which would indicate that the low expression level obtained for some of them (#11 and #9) might be related to poor endocytosis.







# Case Study: LNP formulation screening for Cancer Cells Delivery

## Conclusion

- We managed to select from 15 different formulations with different lipids compositions for their efficiency *in vitro* in HCT116 cell lines.
- Adding fluorescent lipids showed that some formulations fused with cellular membranes but did not lead to protein translation, suggesting some defects during endocytosis.
- This screening allowed us to identify the most potent formulations that would be used for transfection of cancer cells (HCT116) cells.


	4	5	12
EE (%)	63,3	99,6	84,6
z-average (nm)	74,8	120,1	104,3
PI	0,14	0,23	0,07

BF			
eGFP			

## Added Value

- Access to different innovative lipids (ionizable, and others) and to provide ways of tracking LNPs.
- We offered an integrative approach, from LNP preparation to downstream analyses.
- These formulations constitute a cornerstone for *in vivo* experiments.



# mRNA delivery – Our services


Encapsulate higher mRNA quantities



## Upscaling

- We used **off-the-shelf lipid mixtures** (LNP-0315, ABP Biosciences) **and eGFP mRNA** (Cap1, 120A, NI-Methyl-Pseudouridine) **for simplicity purpose.**
- We tested different formulation parameters.

	Flow rate ratios			N/P ratios			Total Flow rate (ml/min)		
	3	4	5	4	6	10	1	2	4
z-average (nm)	156,8	149,7	168,1	155,8	173,0	149,9	183,9	186,0	163,9
Polydispersity index	0,014	0,024	0,033	0,037	0,027	0,017	0,040	0,047	0,020
% encapsulation	89,07	74,78	82,33	87,17	87,12	70,27	77,52	83,31	87,32




- Transfection efficiency was assessed in different cell types (Fibroblasts, MCF, and HCT116) with similar results.
- The next step is to formulate higher quantities of mRNA using the formulations selected in the previous step.



# Webinar Agenda

1. | Tebubio at a glance
2. | Lipid NanoParticles basics
3. | Delivery & Formulation Challenges
4. | Overview of Tebubio Delivery Platform
5. | Case Study: LNP formulation screening for Cancer Cells Delivery
6. | Innovations through our parnters
7. | Live Q&A Session

# Access to a broad range of validated Lipids



# Access to a broad range of validated Lipids



## Network of high-quality partners

- Access to a broad range of validated lipids
- Stay up to date with last development
- Provide a solution to all your projects

## In vitro validation and targeted delivery



Performance of best composition per LipidBrick® lipid

	IM3c	IM12c	IM13c	IM15c	IM16c	IM21.7c	IM22c	IM25c
Hep G2	Orange	Green	Green	Pink	Green	Orange	Orange	Green
HEK-293	Green	Green	Green	Green	Green	Orange	Green	Green
C2C12	Pink	Green	Green	Orange	Green	Orange	Orange	Green
Huh-7	Green	Green	Green	Green	Green	Green	Green	Green
A-498	Orange	Green	Green	Orange	Orange	Pink	Green	Green
Caco-2	Pink	Orange	Pink	Orange	Orange	Orange	Orange	Green
T cells	Orange	Green	Green	Green	Green	Green	Green	Orange

Extensive library of ionizable/cationic lipids from the literature.


# Access to a broad range of validated Lipids



## Network of high-quality partners

- Access to a broad range of validated lipids
- Stay up to date with last development
- Provide a solution to all your projects

## Ex Vivo Studies



*LipidBrick®IM21.7c-LNPs allows Efficient & Durable CAR Expression in Primary T cells using low mRNA amount (~0.8 µg CAR mRNA/106 T cells) & without the need of ApoE supplementation.*

**LipidBrick®** efficient tool for CAR expression

# Access to a broad range of validated Lipids



## Network of high-quality partners

- Access to a broad range of validated lipids
- Stay up to date with last development
- Provide a solution to all your projects





## mRNA based Vaccines




Off-the-shelf formulation (ALC-0315, SM-102)



Innovative lipids for efficient and safe vaccines solutions



### VIRAL TITER



**Our partners have performed extensive analyses of their lipids efficiency for vaccines *in vivo*.**

# Access to a broad range of validated Lipids



## Network of high-quality partners


- Access to a broad range of validated lipids
- Stay up to date with last development
- Provide a solution to all your projects

## In Vivo Targeting




Polyplus<sup>+</sup>

**LipidBrick<sup>®</sup>** efficient lipid for extra-hepatic targeting, in itself or with other ionizable lipids.



# certest

### In vivo TISSUE SELECTIVITY



All Certest lipids have been validated *in vivo*

# Support you to identify the best solution for your project

## Take home message

### **Propose multiple answer to your projects**

Thanks to our extensive lipids' library : off-the-shelf formulations, "classical" ionisable lipids, innovative ionisable/cationic lipids, conjugated lipids, etc.

### **Flexibility around nucleic acids encapsulation**

Whether you want to find the right formulation or want to validate *in vivo* your targets, we can offer adapted solutions.

### **Finding the right formulation for your projects**

through our established partners, we have access to multiple approaches for *in vivo*, *in vitro*, *ex vivo* or mRNA-based vaccines applications

### **Comprehensive Services**

Our experts will help you to move from proof-of-concept studies to validation of your strategies and find the solution to your problematic.

# Webinar Agenda

1. | Tebubio at a glance
2. | Lipid NanoParticles basics
3. | Delivery & Formulation Challenges
4. | Overview of Tebubio Delivery Platform
5. | Case Study: LNP formulation screening for Cancer Cells Delivery
6. | Innovations through our parnters
7. | Live Q&A Session





tebubio

Facilitators of Life Sciences Research

Any Questions ?





tebubio

Facilitators of Life Sciences Research